

# HISTO-PATHOLOGY SOLUTION EXPERT



**BEST PRODUCTS**

**♥ TOP SERVICE**



# We Provide

Wexis provide Histo-pathology Solution Expert (HSE) assist to resolve and optimize the pathology laboratory management.

HSE- This comprehensive solution is designed to optimize and improve the operational efficiency of routine pathology, including Product supply, distribution, operations, and value-added services.



## Faultless Product Supply Chain

Provide a variety of products, a full set of routine pathological product supply, able to fulfill the experimental requirements.



## Professional Distribution Services

Simplify the supply system, and to achieve a unified distribution, delivery time under control.



## Value-added Services

Provide targeted, personalized value-added services program, significantly reducing the workload of pathology technicians.

+120%  
Productivity Increase

## Customer Advantages

- Maximize the operations efficiency
- Decrease operational risk
- Take advantage of Industry Best Practices
- Provide value-added services



### Purchasing and Inventory Management

One-stop purchasing, reduce costs and workload of the procurement staff. Assist inventory management, and real-time replenishment.



### Standardized implementation

Standardized the Pathology Producing Process; Quality control.



### Cost Accounting

Pathological production cost of regular accounting; the basic realization of zero inventory, reduce capital occupation



# Service Scheme

Lifting worries of technical staff through the provision of mature, thoughtful and meticulous routine services, so that they are able to throw themselves into the pathology of professional and technical work.

## Product Supply

Form Grossing → Mounting, nearly 20 specifications of consumables / reagents will be used. One-stop purchase, fully meet the needs of pathological products

## Quality Assurance

The products under the medical device management are obtaining the registration certificate. Quality assurance, and traceability.



## Procurement and Logistics

Service staff docking, one-stop procurement, according to the actual situation to adjust the distribution frequency, and by hand or car unified distribution.

## Cost and Settlement

Regular check and settle the costs of consumables, bills centralized management to reduce the procurement staff and financial staff workload.

## Butler service

Assist to manage inventory to avoid excess / small inventory of the occurrence of real-time replenishment. Follow up daily management of consumables, laboratory preparation and so on. Data sharing, assist in data extraction and analysis.

## Standardized Operation

Partly / completely replace the manual preparation of reagents to liberate medical workers away from chemical gas; maintain standardized, stabilize the level of pathological specimens and sections of the production, staining and diagnostic.

## Service Advantages



# Solution

No.	Steps	Reagents/ consumables	Model / Specifications	Cost
1	Grossing	Trimming Knife Handle	Feather's Handle	
2		Trimming Knife	Feather's Trimming Knife	
3		Embedding Cassettes	Wexis 2000pcs/ carton	
4	Fixation	Tissue Fixative	Wexis 10%/ 13% Neutral Buffered Formalin, 15ml-20L provided	
5	Dehydration	Ethanol (75%/85%)	500ml/bottle, 20bottles/box	
6		Ethanol (95%/100%)	500ml/bottle, 20bottles/box	
7	Cleaning	Xylene	500ml/bottle, 20bottles/box	
8	Infiltration	Paraffin Wax(54-56°C)	Wexis 8KG/carton	
9		Paraffin Wax (57-60°C)	Wexis 8KG/carton	
10	Embedding	Stainless Steel Embedding Mould	Wexis 4 types of Stainless Steel Embedding Mould	
11		Paraffin Wax	Wexis 8KG/carton	
12	Sectioning	Disposable Blades	Feather R35/N35/A35, 50pcs/box	
13		Microscope Slides	Wexis 50pcs/box, 50boxes/carton	
14	Staining	Concentrated Salt	500ml/bottle	
15		Hematoxylin Solution	Wexis 500ml/bottle, 6bottles/carton	
16		Eosin Solution	Wexis 500ml/bottle, 6bottles/carton	
17		Special Stain Solution	Wexis 27 formulas can be choose.	
18	Mounting	Mounting Medium	100 ml/bottle	
19		Cover Glass	Wexis 24*50mm 100pcs/box	



# Company Profile

**Wexis** Group Limited was established in 2009, and as an enterprise provides professional production and operation of pathological consumables, reagents and equipment, and we always committed to the optimization of pathological processes and customer service improvement. Currently, our products and services have been covered 26 provinces, municipalities and autonomous regions, Hong Kong & Macao S.A.R., South Korea, Thailand, Russia, Sri Lanka, UAE, India, Australia, Netherlands, UK and Austria.

Our products line covered pathological equipments reagents and consumables, including Floatation Drying Workstation, Wax Dispenser, Wax Trimmer, Cold Plate, Tissue Embedding Center, Automatic Slides Stainer, Rotary Microtome, Cryo-Microtome, Air Purifier, Reagent Recycler, Grossing Station and Specimen Filing System; Tissue Fixatives, H&E Solution, Special Stain Solution, Paraffin Wax, Microtome Blades(Feather), Embedding Cassettes, and Microscope Slides etc.

We have established a comprehensive service system, distribution system and storage capacity in Guangdong Province and over China.





# Be Your Pathological Manager

We Follow We Understand We Serve



CONTACT



OFFICE

**Head Office:**  
**Wexis Group Limited**  
**Address: Unit 1611, 16/FI,**  
**Yuen Long Trading Centre,**  
**33 Wang Yip Street, Yuen Long,**  
**NT, Hong Kong**  
**Tel: 852-2567 2882**  
**Fax: 852-2567 6333**

**China Office:**  
**Guangzhou Wexis Biotech Limited**  
**Address: No.4 of Yin Sha Road 2,**  
**Wexis Industrial Park,**  
**Xin Tang Sha Pu, Zeng Cheng**  
**District, Guangzhou, P.R. China**  
**Zip-Code: 511338**  
**Tel: 8620-8779 3718 Fax: 8620-8779**

**Website: [www.wexisgp.com](http://www.wexisgp.com)**  
**Email: [marketing@wexisgp.com](mailto:marketing@wexisgp.com)**

