

Wexis

Wexis-Pro Embedding Cassettes

- P.P Material; Stripe Holes
- Unique design with sliding slot, easy assembly & disassembly of the lids
- Both machine printing & handwritten marking possible
- Quality materials, acid, alkali and corrosion-resistant



*Can be used with our optional Stainless Steel Cassette Lid

| Color | <input type="checkbox"/> White | <input type="checkbox"/> Blue | <input type="checkbox"/> Purple | <input type="checkbox"/> Yellow | <input type="checkbox"/> Green | <input type="checkbox"/> Pink | Packing | Description |
|-----------|--------------------------------|-------------------------------|---------------------------------|---------------------------------|--------------------------------|-------------------------------|---|---------------------------------|
| Order No. | 2200001 | 2200004 | 2200007 | 2200010 | 2200013 | 2200016 | 500 pcs/box; 8 boxes/carton | Embedding Cassette, without lid |
| | 2200003 | 2200006 | 2200009 | 2200012 | 2200015 | 2200018 | 500 pcs/box; 8 boxes/carton Cassettes and lids are packed separately | Embedding Cassette, with lid |

Wexis-Pro Embedding Cassettes in plastic transparent tubes

- Working with cassette labeler, improves work efficiency;
- Quality materials of the tubes, deformation-resistant & smooth transport



| Color | <input type="checkbox"/> White | <input type="checkbox"/> Blue | <input type="checkbox"/> Purple | <input type="checkbox"/> Yellow | <input type="checkbox"/> Green | <input type="checkbox"/> Pink | Packing | Description |
|-----------|--------------------------------|-------------------------------|---------------------------------|---------------------------------|--------------------------------|-------------------------------|--|--|
| Order No. | 2200022 | 2200023 | 2200024 | 2200025 | 2200026 | 2200027 | 75 pcs/tube; 60 tubes/carton | Transparent Tubed Cassettes, without lid |
| | 2200028 | 2200029 | 2200030 | 2200031 | 2200032 | 2200033 | 75 pcs/tube; 60 tubes/carton Cassettes and lids are packed separately | Transparent Tubed Cassettes, with lid |


Stainless Steel Embedding Mould

- Quality stainless steel, durable;
- Four different specifications, deepening design in sample trough, suitable for all types of specimens of different sizes

| | | | | | |
|-----------|-------------------------|-----------|-----------|-----------|--------------------------|
| Order No. | 8202101 | 8202201 | 8202301 | 8202401 | 8202001 |
| Size | 37*24*6mm | 30*24*6mm | 24*24*6mm | 10*10*5mm | Set of 4 sizes |
| Packing | 10pcs/pack, 10packs/box | | | | 10sets/pack, 10packs/box |




Standard Embedding Cassettes



- Stripe Holes, with removable lid
- P.O.M. Material; can be used on cassette printing machines

| | | | | | | | |
|-----------|-----------------------------|----------|----------|----------|----------|----------|----------|
| Order No. | 2200037W | 2200037Y | 2200037G | 2200037P | 2200037B | 2200037T | 2200037M |
| Color | White | Yellow | Green | Pink | Blue | Tan | Purple |
| Packing | 250 pcs/box; 8 boxes/carton | | | | | | |



- Round Holes, without lid
- P.O.M. Material; can be used on cassette printing machines

*Can be used with our optional Stainless Steel Cassette Lid

| | | | | | | | |
|-----------|-----------------------------|----------|----------|----------|----------|----------|----------|
| Order No. | 2200038W | 2200038Y | 2200038G | 2200038P | 2200038B | 2200038T | 2200038M |
| Color | White | Yellow | Green | Pink | Blue | Tan | Purple |
| Packing | 500 pcs/box; 8 boxes/carton | | | | | | |



- Stripe Holes, with detachable lid, Assembled
- P.O.M. Material; can be used on cassette printing machines

| | | | | | | | |
|-----------|--|----------|----------|----------|----------|----------|----------|
| Order No. | 2200039W | 2200039Y | 2200039G | 2200039P | 2200039B | 2200039T | 2200039M |
| Color | White | Yellow | Green | Pink | Blue | Tan | Purple |
| Packing | 500 pcs/box, 4 boxes/carton; Cassettes and lids are packed separately | | | | | | |



- Stripe Holes, without lid
- P.O.M. Material; can be used on cassette printing machines

*Can be used with our optional Stainless Steel Cassette Lid

| | | | | | | | |
|-----------|-----------------------------|----------|----------|----------|----------|----------|----------|
| Order No. | 2200040W | 2200040Y | 2200040G | 2200040P | 2200040B | 2200040T | 2200040M |
| Color | White | Yellow | Green | Pink | Blue | Tan | Purple |
| Packing | 500 pcs/box; 8 boxes/carton | | | | | | |



Stainless Steel Cassette Lid

- Quality Stainless steel, reusable, economical & environmentally friendly;
- Compact & easy loading and disassembly

| | |
|-----------|------------|
| Order No. | 2200035 |
| Packing | 100pcs/Box |



Wexis Group Limited
Address: Unit1611, 16/FI, Yuen Long Trading Centre,
33 Wang Yip Street, Yuen Long, NT, Hong Kong
Tel: 852-2567 2882 Fax: 852-2567 6333

Guangzhou Wexis Biotech Limited
Address: No.4 of Yin Sha Road 2, Sha Pu Yin Sha Industrial Park,
Xin Tang, Zeng Cheng District, Guangzhou, P. R, China
Zip-Code: 511388 Tel: 8620-8779 3718 Fax: 8620-8779 3728
www.wexisgp.com



WX-CY-H011401(ENG)
*Image is for illustration purposes only.