



**9L FRESHTAKE™
(2.5 Gallons)**

**OUR BRAND PROMISE:
DOING GOOD WHILE SAVING MONEY™**

**We Make it Easy to Stop
Pouring Your Waste and
Money Down the Drain.
Recycle, Reuse and Realize
Sustainable \$avings!**

A SOLVENT RECYCLER FOR XYLENE, XYLENE-SUBSTITUTES, ALOCHOL & ACETONE

Pure \$avings

- Virtual elimination of solvent purchasing and disposal costs
- 90% recovery from typical histology waste

Ease of Use/ Convenience

- Fully-automated cycles via upgraded touch screen technology
- Smaller footprint dimensions and ergonomic design
- Easier access to a reusable waste bag in a removable spill tray
- No water connection or special ventilation required
- No changing lines or pots
- No cool-down required between runs
- Automatic waste drain and decontamination processes

Safety

- Built-in double spill containment
- No vapor released due to closed-loop ventilation
- ETL listed to UL61010 and UL2208 and CSA

Sustainability

- Allows reuse of solvents for the long-term
- Recovered solvent purity is greater than 99%
- Limied operator labor with automated flush and clean cycles



**9L FRESHTAKE™
Model: MSLV-14-CE-9-A**



**9L FRESHTAKE™
(2.5 Gallons)**

Less Space Required to Recycle in your Laboratory

Footprint Dimensions
42cm(Wide)* 92cm(High)* 61cm(Deep)

Quick Cycle Times

Xylene, Xylene Substitutes & Acetone = 15-20 min per liter/ 3-4 liter per hour
Alcohol=20-30minutes per liter / 2-3 liter per hour

19L Model Also Available