

# WEXIS MICROSCOPE GLASS SLIDE/ COVER GLASS

# Wexis

Specialized manufacturing with internationally recognized standards.


Premium grade thin glass sheet material, precise cutting and unique refining process, guaranteed even thickness and pre-cleaned for every slide. Offers the best optical image under the microscope with precise sizes ideal for use with automatic instruments.

- Fine polished grounded edges, 45/90 degree corners, ease and safe handling
- Even thickness of 1.0-1.2mm, offers users to avoid unnecessary refocusing during microscope observations.
- Variety of models for your choices

## WEXIS PRO Adhesion Microscope Slides


- Super White Glass Class material, pre-cleaned, specialized manufactured for IHC and other cytological and histological tests
- Glass slides surface with special process to offer permanent and more density positive charges to achieve higher adhesion, ideal for securely attaching histological sections and cytological samples onto the slide.
- User choices of M-HPB or M-HPL adhesive properties slides.

Order No.	Description	Packing	Color
62018105PA	Adhesion Microscope Slides Super White Glass Material; Positive Charged Single Painted;	50Pcs/Box; 20Boxes/Carton	White
62018105PAL	Ground Edges; 45° Corners (Optional 90° Corners)	72Pcs/Box; 20Boxes/Carton	



## WEXIS PRO Microscope Glass Slides

- Super White Glass Class material, pre-cleaned, absolutely free of background color, ideal for microscope observation
- 5 painted colors offer ease of categorization and identifications of specimens
- Enameled color coated, provides ideal surface for ease of hand-writing and clearer auto-printer marking for permanent records



Order No.	Description	Packing	Color
62017109PW 62017109PG 62017109PR 62017109PB 62017109PY	Microscope Glass Slide - Processional Grade; Super White Glass Material; Single Painted; Ground Edges; 45° Corners (Optional 90° Corners)	50Pcs/Box; 20Boxes/Carton	White Green Pink Blue Yellow
62017109PWL 62017109PGL 62017109PRL 62017109PBL 62017109PYL		72Pcs/Box; 20Boxes/Carton	White Green Pink Blue Yellow

\*\*Image is for illustration purposes only.

## WEXIS EXCEL Microscope Glass Slides

- Super grade green glass material, pre-cleaned, single/double frosted, suitable for hand writing and printer marking

Order No.	Description	Packing
62017105SP	Super grade green glass material Single Frosted; Ground Edges; 45° Corners(Optional 90° Corners)	50Pcs/Box; 20Boxes/Carton
62017105SPL		72Pcs/Box; 20Boxes/Carton
62017107SP	Super grade green glass material Double Frosted; Ground Edges; 45° Corners(Optional 90° Corners)	50Pcs/Box; 20Boxes/Carton
62017107SPL		72Pcs/Box; 20Boxes/Carton



## Cover Glass

### Borosilicate Glass 3.3 Material

- The premium quality of the borosilicate glass 3.3 material plus the highly sophisticated manufacturing process efficiently avoid the possible attacked by humidity, mold erosions and cover glass sticking problems. Especially suitable for use in most automatic coverslipping machines.

Order No.	Description	Size	Packing
62021200P	Borosilicate Glass 3.3 Material Thickness:0.13-0.16mm	24*24mm	200Pcs/Box; 200Boxes/Carton
62022100P		24*32mm	100Pcs/Box; 200Boxes/Carton
62023100P		24*40mm	
62024100P		24*50mm	
62025100P		24*60mm	



### Super White Glass Material

- Unique premium super white glass material, highly clear and transparent.

Order No.	Description	Size	Packing
62021200	Super White Glass Material Thickness:0.13-0.16mm	24*24mm	200Pcs/Box; 200Boxes/Carton
62022100		24*32mm	100Pcs/Box; 200Boxes/Carton
62023100		24*40mm	
62024100		24*50mm	
62025100		24*60mm	

Wexis Group Limited

Address: Unit1611, 16/FI, Yuen Long Trading Centre,  
33 Wang Yip Street, Yuen Long, NT, Hong Kong  
Tel: 852-2567 2882 Tel: 852-2567 6333

Guangzhou Wexis Biotech Limited

Address :No.4 of Yin Sha Road 2, Sha Pu Yin Sha Industrial Park,  
Xin Tang, Zeng Cheng District, Guangzhou, P. R, China  
Zip-Code: 511388 Tel: 8620-8779 3718 Fax: 8620-8779 3728