

# Wexis

## DR230 Slide Drying Hotplate

Large  
hotplate surface

Slim and light  
Design



Wexis DR230 Slide Drying Hotplate is a key unit offering glass slide drying function for histology. It is always used together with the floatation bath.

\* Utility Model Patent Number: ZL 2012 2 0591789.9  
Design Patent Number: ZL 2012 3 0519819.0

### Built-in Micro-processor Controlled Heating System

- Maintains the required slide drying temperature set by the operator, effectively dries the water on the glass slides after floatation, improving quality of slide staining.

### Unique "Stacking Tiles" design

- The 4 beveled sections hotplate enables higher drying efficiency and enough room for users taking up or placing the glass slides easily onto the hotplate. Operator is thus protected from accidental touching the heated hotplate.

"Stacking Tiles" Style  
4 beveled sections  
hotplate

Micro-processor  
Controlled  
Heating System

Wexis Group Limited is a global pathology supplies and medical equipment manufacturer. We have advanced professional equipment, dedicated research development laboratory, and production workshop, and are one of the leading manufacturers of pathology supplies in China.

Drying hotplate with large capacity dries a maximum of 48 slides at one time.



Special gill type design PE outer casing avoids water splashes and unit is ideally insulated, offers safe and relaxed operation



LED digital display, temperature settings from 00.0 to 99.0C with auto memory of last settings



Accurate temperature controls with automatic adjustments according to set temperatures

## Technical Data

Dimensions	415mm*310mm*120mm
Dryer Hot Plate Dimension	295mm*270mm
Weight	6.2Kg
Line-in(Voltage)	AC220~50Hz 0.3A
Power Rating	300W
Temperature Display Range	0°C - 99°C
Default Temp. Values of	H: 70°C ; L: 68°C
Dryer Temp. Set	Room Temp - 99°C
Accuracy	Set temp. $\pm 2^\circ\text{C}$
Duration of Reach Set Time	10-15mins
Specifications take at	25°C Room Temp.
Certification	CE Quality Certified Organization



## Ordering Number

Order Number	230-00-00
Model	DR230
Name	Slide Drying Hotplate DR230
Includes	Main Unit Power Cable (1) Instruction (1)
Remarks	Photos in this catalogue are for illustration only.

## WEXIS SMALL INSTRUMENT SERIES



Wax Trimmer TR230



Floatation/Drying Workstation FL220



Wax Dispenser WD360

Manufacturer: WEXIS GROUP LIMITED  
Address / Unit 1611, 16/F, Yuen Long Trading Centre, 33 Wang Yip Street, Yuen Long, N.T, Hong Kong  
Tel / 852-2567-2882 Fax / 852-2567-6333

China Head Office: GUANGZHOU WEXIS BIOTECH LIMITED  
Address / No.4 of Yin Sha Road 2, Sha Pu Yin Sha Industrial Park, Xin Tang, Zeng Cheng District,  
Guangzhou, P.R, China (Posted Code: 511388)  
Tel / 8620-8779 3718 Fax / 8620-8779 3728

[www.wexisgp.com](http://www.wexisgp.com)