

Wexis



Cryo-Cabinet

- Special design for pathological biopsy specimen storage
- Dual compressor refrigeration system, double protection, prevents formaldehyde volatilization effectively
- Vacuum-heated glass door, no frosting
- Double open door design with big capacity storage space, mobile specimen storage plate, easy-to-operate

Order No.:	736-00-00
Model No.:	SF736
Product Name:	Cryo-Cabinet SF736
Dimensions (LxWxH):	1200mm*650mm*1910mm
(Line-in) Voltage:	AC220V~50Hz
Power Rating:	240W
Temperature Control:	0 ~ 12°C (Adjustable)
Capacity:	1200L, 20 P.P storage plates for specimen use



Wexis Group Limited

Address: Unit 1611, 16/FI, Yuen Long Trading Centre,
33 Wang Yip Street, Yuen Long, NT, Hong Kong

Tel: 852-2567 2882 Fax: 852-2567 6333

Website: www.wexisgp.com

Email: sales@wexisgp.com