

Wexis Rotary Microtome Series




Safe, convenient, comfortable operation function

Advanced and precise micro system for specimen advance

Ergonomic and unique design

Superior and safe sectioning function



Wexis Rotary Microtome Series



Rotary Microtome M3055

Order No.:	3055-00-00
Model No.:	M3055
Dimensions (LxWxH):	560*470*320mm
Net Weight:	35kg
Max Specimen Size:	60*50*40mm
Coarse feed:	Manual operation through coarse handwheel
Specimen Positioning:	Horizontal 8°, vertical 8°
Specimen Retraction:	About 60µm
Section Thickness Range:	1~60µm (best sectioning thickness 3~5µm)
Selected Thickness:	1µm~10µm, 1µm increments 10µm~20µm, 2µm increments 20µm~60µm, 5µm increments Trim Thickness: 10µm and 30µm



Rotary Microtome M3050

Order No.:	3050-00-00
Model No.:	M3050
Dimensions (LxWxH):	550*465*280mm
Net Weight:	25kg
Specimen Horizontal Stroke:	23mm
Specimen Vertical Stroke:	53mm
Specimen Clamp Adjustable Angle:	8°
Blade Horizontal Stroke:	25mm
Section Thickness Range:	1~60µm
Selected Thickness:	0.5µm~2µm, 0.5µm increments 2µm~10µm, 1µm increments 10µm~20µm, 2µm increments 20µm~60µm, 5µm increments Thickness Deviation: ±20%



Wexis Group Limited

Address: Unit1611, 16/FI, Yuen Long Trading Centre,
33 Wang Yip Street, Yuen Long, NT, Hong Kong
Tel: 852-2567 2882 Fax: 852-2567 6333

Guangzhou Wexis Biotech Limited

Address: No.4 of Yin Sha Road 2, Sha Pu Yin Sha Industrial Park,
Xin Tang, Zeng Cheng District, Guangzhou, P. R, China
Zip-Code: 511388
Tel: 8620-8779 3718 Fax: 8620-8779 3728

Website: www.wexisgp.com E-Mail: sales@wexisgp.com

*Image is for illustration purposes only