



DANGALARK
Pathology & Histology - Supplies & Technology


***The modular
space saving
cost effective
solution
to
cassette & slide
storage***

DANGALARK
STORAGE SYSTEM

With the capacity of 1200+ slides or 420+ cassettes per drawer, you've found the solution to your laboratory's archiving & storage needs.



- Endless stacking potential - *from small to large as needed!*
- Colour code your storage with five colours of drawers - *separate your lab sites or work year.*
- Mobile metal cabinet available for easy transport of drawers from workplace to archives.
- High visibility drawer face allows both non-stick and sticky labels - *identify your storage easily.*
- Pangalark Drawer & Sleeve Dimensions - 12cm W x 56cm L x 10.5cm H



THE LARK DRAWER

Standard slides and cassettes are stored in rows of three.
1,200+ standard glass slides or
210+ standard cassettes
capacity per drawer.



THE LARK INSERT

The Lark Insert Drawer, allows you to double the cassette storage of a Lark Drawer to 420+ standard cassettes.



THE LARK SLEEVE

The Sleeve houses the drawer and connects to other sleeves to create the modular stack. Base version shown above is required to start a verticle stack.

Australian Models shown above. International version will have fixed dividers & label holder.

Distributor

Wexis Group Limited
Address: Unit1611, 16/FI, Yuen Long Trading Centre,
33 Wang Yip Street, Yuen Long, NT, Hong Kong
Tel: 852-2567 2882
Fax: 852-2567 6333
Website: wexisgp.com
E-Mail: sales@wexisgp.com

Wexis