

WEXIS TR302 Wax Trimmer

For Paraffin Block Trimming



WEXIS™ Trimmer TR302

—————Your partner
to offer accurate
paraffin block trimming and repairing
for high quality microtomy.

*Patent Number: ZL201220098285.3

Wide top heating plate

- Wide plain area, allows initial coarse melting of excess wax around specimen blocks.
- Allows possible trimming of 5-6 blocks at one time

Slanted Plate Surface and Wax Scrap Panel

- With sufficient and consistent heating for efficient melting of excess wax on blocks.
- With shallow raked groove design to provide smooth trimming and efficient drainage of excess molten wax

“Drawer Type” Removable Wax Collection Box

- With replaceable liner, waste wax is efficiently collected and easily disposed of.
- Workbench is kept clean and clear of splashed waste wax.



WEXIS Group Limited is a global pathology supplies and medical equipment manufacturer. We have advanced professional equipment, dedicated research development laboratory, and production workshop, and are one of the leading manufacturers of pathology supplies in China.

Wexis

PRODUCT FEATURES

Description & Order Number

Dimensions	285x160x120 (mm)
Weight	1.95kg
Line-in (Voltage)	AC230V 50Hz
Rating	50W
Surface Temp.	78°C (±2°C)
Heat-up time	5-8 mins
Working Temp.	5°C - 40°C
Wax collection drawer capacity	300ml
Order No.	302-0000
Model No.	TR302
Description	Wax Trimmer TR302
Includes	Main Unit (1) Waste Wax Collection Box (1) Wax Disposal Trough (10pcs) Power Cable (1) Instruction (1)

Certification **CE**
Quality Certified Organization
CTS



Low profile; low center of gravity
• keep stability during trimming

90 degree side edge
• offers trimming of small specimen blocks possible

Grip-type hand-rest design
• offers stable and easy control during multiple cassettes trimming

Special plastic body
• waterproof, insulated and heatproof



WEXIS SMALL INSTRUMENT SERIES



Floatation/Drying Workstation FL220



Slide Drying Hotplate DR230



Wax Dispenser WD360

Manufacturer: WEXIS GROUP LIMITED
Address / Unit 1611, 16/F, Yuen Long Trading Centre, 33 Wang Yip Street, Yuen Long, N.T, Hong Kong
Tel / 852-2567-2882 Fax / 852-2567-6333

China Head Office: Guangzhou Wexis Biotech Limited
Address / No.4 of Yin Sha Road 2, Sha Pu Yin Sha Industrial Park, Xin Tang, Zeng Cheng District,
Guangzhou, P.R, China (Zip-Code: 511388)
Tel / 8620-8779 3718 Fax / 8620-8779 3728

