

Wexis


Pathology Grade

Paraffin Wax

- Good toughness, offers ultra thin and smooth section with no wrinkles
- Easy dewaxing , ideal staining result
- High purity, with no odor



Order Number	2200045	2200046
Melting Point	54~56°C	57~60°C
Oil Content	Max 0.5%	Max 0.5%
Pin Penetration Depth(25°C)	16~20	16~20
Packing	1KG/pack, 8 packs/carton	1KG/pack, 8 packs/carton
Color	White	White



Wexis Group Limited

Address: Unit 1611, 16/F, Yuen Long Trading Centre,
33 Wang Yip Street, Yuen Long, NT, Hong Kong

Tel: 852-2567 2882

Fax: 852-2567 6333

E-mail: sales@wexisgp.com